





Analytical Solutions

Filter Vials, Empty Columns, **Collection & Filter Plates**


2091 Rutherford Rd., Suite 150

- Carlsbad, CA 92008
 - info@htslabs.com 🗖
 - 800 541.4792 🕓
 - 760 757.8080 📞
 - htslabs.com 🔇

TIC-PL-082-201 Rev. B



htslabs.com



Analytical Solutions

Filter Vials & Empty Columns

	A Little About Thomson
1	Solutions At Work [™]
2	Open to Collaboration
	About Filter Vials
3	An Introduction to Filter vials
4	How Filter Vials Work
5	Filter Vial Membrane
6	What Applications can the Filter Vial be Used For?
7	What do Filter Vials Replace in the Lab?
8	Optimize Your SPE or QuEChERS Workflow
9	A Comparison of the Filter Vial Types
11	Plasticizers Content in Filter Vials Compared to Syringe Filters
12	Increase Signal-to-noise Ratio with eXtremelFV® for More Targeted & Accurate Peaks
13	Chemical Compatibility
14	Compound Compatibility
15	High Viscosity Presses

 Empty Columns

 17
 SINGLE StEP® Empty Columns


	Collection & Filter Plates
19	Well Plates & Filter Plates
21	Screening Protocols for Mammalian/Insect Cells
22	Screening Protocols for E. coli and other Microbes
23	Clarification Protocol for Mammalian/Insect Cells
24	Clarification Protocol for E. coli and other Microbes
	Part Numbers
25	StandardlFilter Vials
25	eXtremelFV®
25	Low EvaplFilter Vials
26	nanolFilter Vials®
26	High Viscosity Presses
27	SINGLE StEP [®] Empty Columns

- **28** Collection & Filter Plates
- 29 Plate Seals/Capmats

A Little About Thomson

SOLUTIONS AT WORK[™]

Thomson sells innovative single-use Solutions At Work[™], our mission is to provide technical expertise while partnering with our customers to deliver practical scientific innovations enabling scientific advancements in pharmaceutical, biotech, environmental/food, toxicology/forensics, and contract manufacturing industries.



Open to Collaboration

INNOVATIVE PRODUCT LINE

Scientists around the world are discovering new ways to use Thomson Filter Vials. Whether testing pharmaceuticals, performing toxicology, or testing for drugs of abuse Thomson Filter Vials have proven to be indispensable tools for sample prep when using HPLC, GC, LC-MS, or GC-MS, methodologies.

Thomson offers a full line of shake flasks and accessories with aboveaverage yields and higher working volumes, designed specifically for insect/mammalian, or microbial/*E. coli* cells based on an understanding and experience of lab operations.

Our well-plate catalog continues to grow and provide the highest quality plates, ready for robotics, cell culture, synthesis, or analysis.

SINGLE StEP $\ensuremath{^{\tiny \circledast}}$ Empty Columns are ready for the addition of sorbents or resins depending on the application.

If you have unique needs or need a new product please reach out to us. We look forward to collaborating with you.



An Introduction to Filter Vials

Thomson Filter Vials are a single system which replaces HPLC Vials, HPLC Caps, Syringes, & Syringe Filters for the filtration of samples. In 15 seconds, Thomson Filter Vials filter samples in an autosampler-ready vial.

Key Features

- Same Size as a standard HPLC Vial and will fit easily into any standard HPLC vial machine or tray
- PTFE, PVDF, PES and Nylon membranes are available depending on the percentage of organic solvent in the sample and the amount of protein binding
- Pore sizes of either 0.2µm or 0.45µm will provide the perfect degree of filtration needed from viscous to clarified samples
- Versatility is built into Thomson's line of Filter Vials. Whether your samples are low volume or viscous or particulate-laden or contain a high volatility organic solvent Thomson has a Filter Vial to fit your needs



Syringe Filter Built In

Equivalent to A Syringe Filter Built Into Your HPLC Vial

Filter Vials are equivalent to a syringe filter built into your HPLC vial. Even samples that appear clear to the eye potentially have particulates that can clog the machine, causing down time and costly maintenance. Filter Vials increase productivity by eliminating a transfer step required when using a syringe filter.

How Filter Vials Work

Similar to How A French Press Works...



Easy As 1, 2, ... Done!

In Two Steps

- 1. Deposit 450µL of sample into shell vial
- 2. Insert plunger into the outer shell & press

15 Seconds

48) can be used.

syringe packaging and add a syringe filter.

Filter Vial Membrane

What Applications Can the Filter Vial be Used For?

Membrane Pore Size

The recommended membrane pore size for sample filtration is based on the cell or cell debris content of the sample and the particle size of the packing material in the chromatography column used to analyze the sample. If the sample contains cells or cellular debris, then a 0.2µm pore size membrane is recommended to maintain system sterility.

Which to use?

- 0.2um Pore Size
- Cells or Cell Debris in Sample
- Chromatography Column Particle Size <3µm
- 0.45µm Pore Size
- Chromatography Column Particle Size >3µm

With Thomson's family of Filter Vials and membranes available to you, finding ways to replace cumbersome and expensive syringe filters in the lab is easy. Here are just some of the documented applications you can use Filter Vials for in your lab today. See our Technical Library at htslabs.com to see a full list of applications. We work hard with small and large companies to produce proven protocols and methods for our products. If you find a use for Filter Vials in your workflow we would love to hear about it.

Membrane Material

The recommended membrane for sample filtration is based on the percentage of organic solvent in the sample and the amount of protein binding.

Compatibility

For chemical or compound compatibility with our Filter Vials & membranes see the Chemical Compatibility Index & Compound Compatibility Index in our Technical Library.

	Aqueous	>50% Organic	Low Protein Binding
PTFE			
PVDF			
Nylon			
PES			

Thomson's Technical Library

You can find application notes, videos and more information on our products by visiting our website at htslabs.com.



	nano Filter Vial®	StandardlFilter Vial	Low EvaplFilter Vial	eXtremelFV®
10µL-250µL				
450µL				
UPLC Compatible				
GCMS Compatible				
30% Particulates				
Viscous				
Replacement for SPE				
General Liquids < 10% particulates				
Cell Fermentation				
Particulate Removal				
Automation Compatible				
Small Molecules				
Food & Supplements				
Toxicology				
Pesticides				
Environmental				







NO MORE Syringes

Indundandandand

NO MORE Syringe Filters

What do Filter Vials Replace in the Lab?

Thomson Filter Vials simplify general filtration by replacing syringes & syringe filters, microcentrifuge spin columns, and/or liquidliquid extractions.

Applications for Thomson Filter Vials include all sample types to be analyzed by HPLC, UHPLC, LC-MS, and GC-MS.




NO MORE HPLC Vials & Caps

A Comparison of the Filter Vial Types

Filter Vial

Standard For Most Samples

Max Fill Vol. 450µL Dead Vol. 120µL




Key Features

- General purpose filtration • <10% particulates
- Pre-slit septum

nan•|Filter Vial.

When Every µL Counts

Max Fill Vol. 250μL Min Fill Vol. 10μL (for 2μL injection)



Key Features

10μL sample for 2μL injection
Available with pre-slit or non-slit septum

Replaces in the lab

- Syringe Filters
- Syringes
- HPLC Vials/Caps

Applications

- •120µL-450µL
- General Liquids < 10% particulates
- Particulate Removal
- Automation Compatible
- Small Molecules
- Food & Supplements
- Toxicology
- Environmental

Replaces in the lab

- Centrifugation & Spin Filters
- Small Volume Syringe Filters
- Syringes
- High Recovery Vials/Caps
 Inserts with HPLC Vials/Caps


Applications

- •10µL-250µL
- General Liquids < 10% particulates
- Cell Fermentation
- Particulate Removal
- Automation Compatible
- Small Molecules
- Toxicology
- Pesticides
- Environmental

EXTREME/FV.

Multi-Layered Filtration

Max Fill Vol. **450µL** Dead Vol. **120µL**



Key Features

Used for Particulate Laden Samples
Contains a Depth Pre-Filter
Pre-slit septum

Replaces in the lab

Syringe FiltersSyringesHPLC Vials/Caps


Applications

- •120µL-450µL
- \leq 30% Particulates
- Viscous
- Replacement for SPE
- Cell Fermentation
- Particulate Removal
- Automation Compatible
- Small Molecules
- Food & Supplements
- Toxicology
- Pesticides
 Environmental

Low Evap | Filter Vial

Standard For Most Samples

Max Fill Vol. 450µL Dead Vol. 120µL



Key Features

- General purpose filtration
- Non-split septum
- <10% particulates
- Evaporation rate <0.4% over 24-hour

Replaces in the lab

- Syringe Filters
- Syringes
- HPLC Vials/Caps

Applications

- •120µL-450µL
- General Liquids < 10% particulates
- Particulate Removal
- Automation Compatible
- Small Molecules
- Food & Supplements
- Toxicology
- Environmental

Plasticizers content in Filter Vials Compared to Syringe Filters

Testing by Takeda Pharmaceutical Company Limited® UPLC - ELSD

Introduction

Thomson Filter Vials are manufactured without the use of plasticizers or mold release agents, making them LC/MS clean. Testing with ELSD, PDA, and MS detection by Takeda Pharmaceutical showed no leaching from Thomson Standard Filter Vial with a 0.45µm, PTFE membrane compared to significant leaching from Millipore Millex-FH® Filter, 0.45µm, hydrophobic PTFE, 4mm. Method: A. Water B. ACN 45-90% with 0.05% TFA Ballistic Gradient over 1.4 minutes using Waters® Acquity[®] UPLC Thomson Filter Vial (patented) Part # 34440 Filter Vial 0.45µm hydrophobic PTFE, w/ Pre-Slit Cap Millipore Syringe Filter Part #:SLFHR04NL Millex-FH® Filter, 0.45µm, hydrophobic PTFE, 4mm, nonsterile.

Method:

A. Water B. ACN 45-90% with 0.05% TFA

Ballistic Gradient over 1.4 minutes using Waters® Acquity® UPLC

Thomson Standard Filter Vial


0.45µm hydrophobic PTFE, w/ Pre-Slit Cap Part#: 34440


Millipore Syringe Filter


Millex-FH® Filter, 0.45µm, hydrophobic PTFE, 4mm, non-sterile. Part #: SLFHR04NL


Plasticizers



















Thomson is not affiliated with Takeda Pharmaceutical Company®, Millipore®, Waters® or their products

Increase Signal-to-Noise Ratio with eXtreme|**FV**[®] for More Targeted & **Accurate Peaks**


Matrix Effects & Ion Suppression:

Analytes are obscured by the matrix like the octopus in this photo is difficult to find among its surroundings.



Low Signal to Noise Ratio

Difficult to find analyte in the matrix



Octopus images courtesy Jukin Video




Strong Signal; Noise Lessened:

By adding compounds to the eXtremelFV® the signal to noise ratio is increased allowing you to find the analyte with ease.



High Signal to Noise Ratio

In this example the addition of C-18 to eXtremelFV® with your sample binds excess compounds to C-18 and the Matrix clears up allowing you to see analyte peaks



Chemical Compatibility

	Housing Materials		Filter Membrane		
	Polypropylene	PTFE	PVDF	PES	NYLON
Acetic Acid (glacial) acid, organic	TST	R	R	R	NR
Acetone ketone	R	R	NR	GNR	R
Acetonitrile (ACN) nitrile	R	R	LTD	NR	R
Ammonium Hydroxide caustic	TST	GR	R	NR	TST
Ammonium Sulfate (saturated) salt, aqueous solution	R	GR	NR	ND	R
Amyl Alcohol alcohol	R	R	R	GR	TST
Benzene HC, aromatic	NR	_	_	—	_
Benzyl Alcohol HC aromatic/alcohol	NR	-	_	_	-
Butyl Alcohol alcohol	R	GR	R	GR	R
Chloroform HC, halogenated	NR	-	_	_	-
Cyclohexanone ketone	NR	-	_	_	_
Dimethyl Sulfoxide (DMSO) sulfoxide	R	R	NR	NR	R
Dimethylacetamide amide	R	GR	NR	NR	NR
Dimethylformamide amide	R	GR	NR	ND	R
Ethyl Acetate ester	TST	R	R	GNR	R
Ethyl Alcohol alcohol	R	R	R	GR	TST
Ethylene Glycol glycol	R	R	R	GR	R
Formaldehyde aldehyde	R	R	R	ND	R
Formic Acid, 50% acid, organic	R	GR	R	ND	NR
Glycerine (Glycerol) glycol	R	GR	R	GR	R
Hexane HC, aliphatic	NR	_	_	_	_
Hydrochloric Acid, 1N (HCL) acid, inorganic	GR	R	R	GR	GR
Hydrochloric Acid, 6N (HCL) acid, inorganic	TST	R	TST	GR	TST
sobutyl Alcohol alcohol	R	R	R	GR	TST
sopropyl Acetate ester	TST	R	R	GNR	R
sopropyl Alcohol alcohol	R	R	R	GR	TST
Lactic Acid, 50% acid, organic/alcohol	R	GR	TST	ND	TST
Methyl Acetate ester	TST	R	NR	GNR	R
Methyl Alcohol alcohol	R	R	R	GR	TST
Methylene Chloride HC, halogenated	NR	-	_	—	-
Nitric Acid, 6N acid, inorganic	TST	R	R	R	NR
Nitrobenzene HC, aromatic	NR	-	_	—	-
Pentane I HC, aliphatic	NR	_	_	—	-
Phenol (aqueous solution) phenol	NR	-	_	—	-
Potassium Hydroxide, 3N caustic	R	R	R	ND	R
Silicone Oils silicone	R	GR	R	ND	R
Sodium Carbonate (aqueous solu-tion) salt, aqueous solution	R	R	R	ND	TST
Nater (Brine) salt, aqueous solution	R	R	R	ND	R
Sodium Chloride (aqueous solution) salt, aqueous solution	R	R	R	ND	R
Sodium Dodecyl Sulfate surfactant/detergent	ND	ND	ND	ND	ND
Sodium Hydroxide, 3N caustic	R	R	R	R	R
Sulfuric Acid (concentrated) acid, inorganic	NR	-	-	-	-
Tetrahydrofuran (THF) ether	NR	_	_	_	_
Toluene HC, aromatic	NR	-	-	-	-
TCA (aqueous solution) acid, organic	R	GR	R	ND	TST
Tween® 20 (aqueous solution) surfactant/detergent	ND	R	TST	ND	TST


R = Recommended | GR = Generally Recommended | NR = Not Recommended | GNR = Generally Not Recommended

LTD = Limited Recommendation | TST = Testing Recommended | ND = No Data Presently Available | — = Not Recommended, polypropylene is NR

Compound Compatibility

	5-Fluorouracil
	18F) Fluoromisondazole, Misiomidazole
	Acetylsalicylic acid
	Alprenolol
ŀ	Amiloride
ŀ	Atenolol
	Azathioprine
	Azodicarbonamide
E	Bleomycin Sulfate
	Caffeine
(Cetirizine
	Chlorothiazide
(Chloramphenicol
	Cimetidine
(Ciprofloxacin
	Cyclosporine A
(Cytarabine
	Diclofenac
	Hydrochlorothiazide
	buprofen
j;	sonicotinic acid
k	Ketamine
L	evofloxacin
L	omefloxacin
N	Metoprolol
N	Mitomycin
Ν	Norphazinamide
Ν	Nadolol
Ν	Nicotinic acid
F	Paclitaxel
F	p-Aminobenzoic acid (PABA)
F	o-aminosalicylic acid
	Pefloxacin
F	Pentoxifylline (PTX)
F	Phenytoin
F	Ranitidine
	Rifampicin
5	Sabeluzole
S	Sulfadozine
5	Sulphasalazine
	Sulpiride
	Ferbutaline
T	Fimolol
T	Franexamic acid
T	Friamcinolone Acetonide
Т	Fropicamide





High Viscosity Presses


The Thomson Filter Vial Press enables high solid content and viscous liquids to be easily filtered through vials. Some fermentation cultures that reach 100OD or particulate laden samples may require the toggle press.

Toggle Press

- Press up to 5 autosampler ready Filter Vials
- Allows for consistency and ergonomic concerns
- Small footprint; sits on bench top
- Works with all Thomson Filter Vials

Multi-Use Press

- Presses up to 48 Autosampler Ready Filter Vials at a time
- Works with 48 position block; block fits some autosamplers
- 48 position block can be transferred to a robot for automation
- Easily Automate Filter Vial Pressing
- Works with all Thomson Filter Vials





Empty Columns Fill different sized columns with a variety of sorbents &

resins for purification application.

Easy to Use

Don't be Limited with Your Column Size

At Thomson, we are aware of the need to customize available apparatus to the individual experiments. Our SINGLE StEP® Empty Columns (patented) provide you the opportunity to fill different sized columns with a variety of sorbents and resins for purification applications. The wide range in which we offer these columns means you are not limited by column size. SINGLE StEP® Empty Columns allow for the simple connection to FPLC/LPLC/MPLC systems.

Key Features

- Acceptable for use with Gravity or FPLC/LPLC/MPLC
- Multiple sizes for scales from 10mL-600mL (4g 300g)
- Durable design for pressures up to 200psi
- Top & bottom connections are standard Luer sizing



50mL

10mL

25mL

Collection Plates

Thomson Well Plates in both 24- and 96-Well configurations are ideal for sample preparation or concentration and feature various well and well bottom shapes to suit your analytical needs. To compliment Thomson well plates we also offer various sealing options including capmats, airporous seals, foil seals and plastic lids.

Well Shape – Square or Round to fit your cell type and culture condition requirements
 Well Bottom Shape – Pyramid, Round and V-bottom to fit your applications
 Well Plate Orientation – Fixed for Robotic Liquid-Handling Systems





Filter Plates

ThomsonFilterPlatesinboth24-and96-Wellconfigurations are designed for analytical sample preparation. Depending on your application we may recommend using a positive pressure manifold, centrifugation or a Thomson Vacuum Manifold.

• Versatility – solid phase extraction and affinity phase adsorption applications involving high throughput robotic Liquid Handling Systems

• Solvent Compatibility – PVDF and PTFE Filter Plates are similar in principle to Thomson Filter Vials but in a 96-well plate

• Long funnel design – Eliminates cross-contamination between sample collection wells by fully inserting below the top of the collection plate

Materials:

• 96-Well Plate, 2mL, Square Well, Pyramid Bottom, Individually Wrapped with Lid Sterile: p/n 931137



Materials:

- 24-Well Plate, 10.4mL, Square Well, Round Bottom, Individually Wrapped with Lid Sterile: p/n 931568
- 24-Well Plate, 10.8mL, Square Well, Pyramid Bottom, Individually Wrapped with Lid Sterile: p/n 931571

96-Well Screening Protocol for Mammalian/Insect Cells

Methods:

- 1. Maintain cell stocks in appropriate growth medium. Split cultures the day before transfection to an appropriate density to ensure log phase growth at the time of transfection.
- 2. Seed cells at 500µL/well. The optimal seeding density will depend on the cell line, please use cell line recommended density.
- 3. Transfect cells according to established transfection protocol. Scale transfection reagent/DNA/feeds on a volumetric basis from what is used for larger scale cultures.
- 4. Seal plates with plastic lids or Airporous Seals and transfer to shaker overnight at 800rpm on a 3mm orbit at 37°C.
- 5. Harvest cells at the time point established for larger scale cultures. Pellet cells by centrifugation at 1000-2500g for 10-20min at 4°C.
- 6. Reserve either the culture media or the pellet depending on the application and proceed to downstream processing.

Notes:

- The most critical factor in cell viability is aeration. Optimal results will be achieved using shakers with 3mm orbit diameters. We do not recommend working in 96well format using shakers with standard 25mm throws.
- Thomson filter plates are a great complimentary product for downstream purification applications.
- •96-Well Filter Plate, 2mL, Long Drip | 25µm Polypropylene: p/n 931919 Maximum suggested centrifugation: 3000g

24-Well Screening Protocol for Mammalian/Insect Cells

Methods:

- 1. Maintain cell stocks in appropriate growth media. Split cultures the day before transfection to an appropriate density to ensure log phase growth at the time of transfection.
- 2. Seed cells at 4-5mL/well. The optimal seeding density will depend on the cell line, please use cell line recommended density.
- 3. Transfect cells according to established transfection protocol. Scale transfection reagent/DNA/feeds on a volumetric basis from what is used for larger scale cultures.
- 4. Cover plates with plastic lids and transfer to shaker overnight at 350rpm on a 12.5mm orbit at 37°C.
- 5. Harvest cells at the time point established for larger scale cultures. Pellet cells by centrifugation at 1000-2500g for 10-20min at 4°C.
- 6. Reserve either the culture media or the pellet depending on the application and proceed to downstream processing.

Notes:

• The most critical factor in cell viability is aeration. For 24-well plates optimal results will be achieved using shakers with 12.5mm. Cultures grown in shakers with standard 25mm throws, will likely need increased rotational speed or decreased culture volume.

96-Well Screening Protocol for E. coli and other Microbes

Methods:

- 1. Pipette 750µL of bacterial growth media containing the appropriate selective antibiotic into each well of a 96-well plate.
- 2. Add the selected colony to each well from either an agar plate or glycerol stock.
- 3. Gently triturate each well manually.
- 4. Seal plates with Airporous Seals and transfer to shaker overnight at 850rpm on a 3mm orbit at 37°C.
- 5. Harvest the plates by centrifugation @ 2500g for 20 minutes.
- 6. Invert the plate to discard the media.
- 7. Process samples according to downstream application (plasmid purification, protein extraction etc.).

Notes:

- •96-well cultures grown in Plasmid+® medium (p/n 446300) typically provide the appropriate biomass for MINI scale plasmid preps.
- The most critical factor in cell viability is aeration. Optimal results will be achieved using shakers with 3mm orbit diameters. We do not recommend working in 96well format using shakers with standard 25mm or 50mm throws.
- If high levels of evaporation are encountered, the well plate & plastic lid (p/n 931134) is recommended to alleviate the issue.
- Thomson Instrument Company's filter plates are a great complimentary product for downstream purification applications. Add appropriate resin.
- •96-Well Filter Plate, 2mL, Long Drip | 25µm Polypropylene: p/n 931919

24-Well Growth Protocol for E. *coli* and Other Microbes

Methods:

- 1. Pipette 4-5mL of bacterial growth media containing the appropriate selective antibiotic into each well of a 24-well plate.
- 2. Add the selected colony to each well from either an agar plate or glycerol stock.
- 3. Gently triturate each well manually.
- 4. Seal plates with Airporous Seals and transfer to shaker overnight at 350-400rpm on a 12.5mm orbit at 37°C.
- 5. Harvest the plates by centrifugation @ 2500g for 20 minutes.
- 6. Invert the plate to discard the media.
- 7. Process samples according to downstream application (plasmid purification, protein extraction etc.).

Notes:

- 24-well cultures grown in Plasmid+® medium (p/n 446300) typically provide the appropriate biomass for MIDI scale plasmid preps.
- The most critical factor in cell viability is aeration. For 24-well plates optimal results will be achieved using shakers with 12.5mm. Cultures grown on shakers with standard 25mm throws will likely need increased rotational speed or decreased culture volume.
- If high levels of evaporation are encountered, use a 24-well plate & plastic lid (p/n 931568) is recommended to alleviate the issue.

Materials:

- 96-Well Plate, 2mL, Square Well, V-Bottom, Raised Lettering | Sterile: p/n 951652C
- 96-Well Plate, 2mL, Square Well, Round Bottom | Irreversible: p/n 931130
- Airporous Seal For Growing Cultures: p/n 899410



Materials:

- 24-Well Plate, 10.4mL, Square Well, Round Bottom, Individually Wrapped | Sterile: p/n 931565-G-1X
- Plus: Airporous Seal for Growing Cultures: p/n 899410
- 24-Well Plate, 10.8mL, Square Well, Pyramid Bottom, Individually Wrapped | Sterile: p/n 931569-G-1X
- Plus: Airporous Seal for Growing Cultures: p/n 899410
- 24-Well Plate, 10.4mL, Square Well, Round Bottom, Individually Wrapped with Lid Sterile: p/n 931568
- 24-Well Plate, 10.8mL, Square Well, Pyramid Bottom, Individually Wrapped with Lid | Sterile: p/n 931571

Materials:

- Thomson Rapid Clear[®] 2, 96-Well 0.2µm Filter Plate, Sterile | CS20 p/n 921746
- Culture plate: 96-Well Plate, 2mL, Square Well, Pyramid Bottom, Individually Wrapped w/ Lid, Sterile | CS20 p/n 931137
- Optional: Airporous Plate Seal, For Growing Cultures, Sterile | Use w/ All Plates | CS100 p/n 899410
- Collection plate: 96-Well Plate, 2mL, Square Well, Round Bottom, Irreversible, Sterile CS20 p/n 931130
- Adhesive Foil Seal | Use w/ All 96- & 24- Well Plates | CS100 p/n 899405-1



96-Well Media Clarification Protocol for Mammalian Cells:

Cell Lines Expi293 and ExpiCHO transient expressing 130kDa, IgG-type protein

Methods:

- 1. Maintain cell stocks in appropriate growth medium with 0.8mL of media / well. Split cultures the day before transfection to an appropriate density to ensure log phase growth at the time of transfection.
- 2. Seed cells at 500µL/well. The optimal seeding density will depend on the cell line, please use cell line recommended density.
- 3. Transfect cells according to established transfection protocol. Scale transfection reagent/DNA/feeds on a volumetric basis as needed.
- 4. Cover plates with plastic lids or airporous seals and transfer to shaker*: A. Expi293 shake speed: 1000 RPM (3mm throw) B. ExpiCHO shake speed: 900 RPM (3mm throw)
- 5. Culture cells for recommended time period.
- 6. Determine viable cell count. Target VCC should be approximately: A. Expi293 7-9 x e6 / mL
 - B. ExpiCHO: 12-15 x e6 / mL
- 7. Pellet cells by centrifugation at 1000-2500 x g for 10-20min at 4°C.
- 8. Assemble Rapid Clear® 2 (RC2) on top of 96-well collection plate.
- 9. Add 400µL of sterile PBS to each RC2 well as rinse.
- 10. Pipette full culture volume, minus cell pellet, into wells of RC2 (700-800µL/ well). Place RC2/collection plate assembly into centrifuge swinging bucket and centrifuge for 5 min at 750-1000 x g.
- 11. Add 600µL of sterile PBS flush to each RC2 well and repeat centrifugation step.
- 12. Disassemble RC2/collection plate assembly and, to prevent evaporation, seal collection plate with foil seal.
- 13. Proceed to Octet column chromatography for protein purification, expect 95% yield (70% yield w/o PBS flush).

Notes:

*The most critical factor in cell viability is aeration. Optimal results will be achieved using shakers with 3mm orbit diameters. We do not recommend working in 96-well format using shakers with standard 25mm or 50mm throws.

96-Well Media Clarification Protocol for Microbial Cells:

E. coli expressing known protein

Methods:

- 1. Maintain cell stocks in appropriate growth medium* with 0.8mL of media/well. Split cultures the day before transfection to an appropriate density to ensure log phase growth at the time of transfection.
- 2. Transfect cells according to established transfection protocol. Scale transfection reagent/DNA/feeds on a volumetric basis as needed.
- 3. Cover plates with plastic lids or airporous seals and transfer to shaker**. Shake at 1000 RPM (3mm throw) for 24 hours.
- 4. Pellet cells by centrifugation at 1000-2500 x g for 10-20min at 4°C.
- 5. Assemble Rapid Clear[®] 2 (RC2) on top of 96-well collection plate.
- 6. Add 400µL of sterile PBS to each RC2 well as rinse.
- 7. Pipette full culture volume, minus cell pellet, into wells of RC2 (700-800µL). Place RC2/collection plate assembly into centrifuge swinging bucket and centrifuge for 5 min at 750-1000 x g.
- 8. Add 600µL of sterile PBS to each RC2 well and repeat centrifugation step. 9. Disassemble RC2/collection plate assembly and, to prevent evaporation, seal collection plate with foil seal.
- 10. Proceed to plasmid purification, expect 95% yield (70% yield w/o PBS flush).

Notes:

*24-well cultures grown in Plasmid+® medium (p/n 446300) typically provide the appropriate biomass for MIDI scale plasmid preps.

**The most critical factor in cell viability is aeration. Optimal results will be achieved using shakers with 3mm orbit diameters. We do not recommend working in 96-well format using shakers with standard 25mm or 50mm throws.

Materials:

- Thomson Rapid Clear[®] 2, 96-Well 0.2µm Filter Plate, Sterile | CS20 p/n 921746
- Culture plate: 96-Well Plate, 2mL, Square Well, Pyramid Bottom, Individually Wrapped w/ Lid, Sterile | CS20 p/n 931137
- Optional: Airporous Plate Seal, For Growing Cultures, Sterile | Use w/ All Plates | CS100 p/n 899410
- Plasmid+[®] Media, Sterile, 1L | CS6 p/n 446300
- Collection plate: 96-Well Plate, 2mL, Square Well, Round Bottom, Irreversible, Sterile CS20 p/n 931130
- Adhesive Foil Seal | Use w/ All 96- & 24- Well Plates | CS100 p/n 899405-1



Part Numbers

Standard|Filter Vial Snap Cap

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm
Cap Color	green	blue	red	yellow	black	pink	grey
Cap Style	snap-cap						
Septum	pre-slit						
Fill Vol.	450µL						
Dead Vol.	120µL						
Part #	35530	35540	35531	35541	35538	35539	35535
Qty/Case	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500

Standard|Filter Vial Screw Cap

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm
Cap Color	green	blue	red	yellow	black	pink	grey
Cap Style	screw cap						
Septum	pre-slit						
Fill Vol.	450µL						
Dead Vol.	120µL						
Part #	34430	34440	34431	34441	34438	34439	34435
Qty/Case	100	100	100	100	100	100	100

eXtreme|FV[®] Snap Cap

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm
Cap Color	green	blue	red	yellow	black	pink	grey
Cap Style	snap-cap						
Septum	pre-slit						
Fill Vol.	450µL						
Dead Vol.	120µL						
Part #	85530	85540	85531	85541	85538	85539	85535
Qty/Case	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500

eXtreme|FV[®] Screw Cap

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm
Cap Color	green	blue	red	yellow	black	pink	grey
Cap Style	screw cap						
Septum	pre-slit						
Fill Vol.	450µL						
Dead Vol.	120µL						
Part #	84430	84440	84431	84441	84438	84439	84435
Qty/Case	100	100	100	100	100	100	100

Low Evap|Filter Vial

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm
Cap Color	green	blue	red	yellow	black	pink	grey
Cap Style	screw cap	screw-cap					
Septum	non-slit						
Fill Vol.	450µL						
Dead Vol.	120µL						
Part #	64430	64440	64431	64441	64438	64439	64435
Qty/Case	100	100	100	100	100	100	100

nano|Filter Vial® Non-Slit

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	NYLON	PTFE	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.2µm	0.2µm	0.2µm
Cap Color	green	blue	red	yellow	black	black	green	grey
Cap Style	screw cap	screw cap	screw cap					
Septum	non-slit	non-slit	non-slit	non-slit	non-slit	non-slit SIL PP	non-slit PTFE SIL PTFE	non-slit
Fill Vol.	250µL	250µL	250µL	250µL	250µL	250µL	250µL	250µL
Dead Vol.	8µL	8µL	8µL	8µL	8µL	8µL	8µL	8µL
Part #	15530	15540	15531	15541	15538	14638	14930	15535
Qty/Case	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	100	100	200 & 500

nano|Filter Vial® Pre-Slit

Membrane	PTFE	PTFE	PVDF	PVDF	NYLON	PES
Pore Size	0.2µm	0.45µm	0.2µm	0.45µm	0.2µm	0.2µm
Cap Color	green	blue	red	yellow	black	grey
Cap Style	screw cap					
Septum	pre-slit	pre-slit	pre-slit	pre-slit	pre-slit	pre-slit
Fill Vol.	250µL	250µL	250µL	250µL	250µL	250µL
Dead Vol.	8µL	8µL	8µL	8µL	8µL	8µL
Part #	25530	25540	25531	25541	25538	25535
Qty/Case	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500	200 & 500

High Viscosity Filter Vial Presses

Press	Description	Capacity	Qty	Part #
Toggle Press	5 Position for Autosampler Ready Filter Vials	5	1	35005
Multi-Use Press	48 Position for Autosampler Ready Filter Vials	48	1	35015

Empty Columns

Description	Case/Qty	Part #
Empty SINGLE StEP® Fritted column w/10 each size: 10mL, 25mL, 50mL, 100mL, 200mL	50	94520-10
Franky CINCLE CEED® Column 10ml or 10 Decompisionith Exit	10	9452086-10
Empty SINGLE StEP® Column 10mL or 4g Reservoir with Frit	100	9452086-100
Empty SINGLE StEP® Column 25mL or 12g Reservoir with Frit.	10	9452088-10
	100	9452088-100
Empty SINGLE StEP® Column 50mL or 25g Reservoir with Frit.	10	9452090-10
	100	9452090-100
Empty SINGLE StEP® Column 100mL or 40g Reservoir with Frit.	10	9452092-10
	100	9452092-100
Empty SINGLE StEP® Column 200mL or 80g Reservoir with Frit.	10	9452094-10
	100	9452094-100
Empty SINGLE StEP® Column 320mL or 160g Reservoir with Frit.	10	9452099
Empty SINGLE StEP™ Column 600mL or 300g Reservoir w/5mL Bottom Resin Reservoir w/Frit	10	9452097-B
Frits		
Head Space Frits for 4g SINGLE StEP® Column	100	491250
Head Space Frits for 12g SINGLE StEP® Column	100	491252
Head Space Frits for 25g SINGLE StEP® Column	100	491253
Head Space Frits for 40g SINGLE StEP® Column	100	491254
Head Space Frits for 80g or 90g SINGLE StEP® Column	100	491256
Head Space Frits for 110g or 160g SINGLE StEP® Column	100	491258
Head Space Frits for 240g or 300g SINGLE StEP® Column	100	491260
Accessories	,	
Dual Ended PP Cap Blue	100	235008

Accessories not sold by Thomson

Tubing Size to Adapt SINGLE StEP® Column to FPLC, AKTA	Upchurch PN#*	Cole-Parmer® PN#*	
1/8th line (obsolete p/n 295821)			
Assembled	LT-215 P-359 P-658 P-655	EW-02022-43 EW-02023-15 EW-02014-14 EW-02014-12	
l/16th line (obsolete p/n 295823)	P-656 P-659	P-656 P-659	

Not all parts may be needed for your set up.

*Thomson Instrument Company is not affiliated with Upchurch, Cole-Parmer® or their products

Vol. Well	Well Shape	Sterility (SAL)	ANSI-SLAS	Ind. Wrap	Compatible with	Capmat/Seal	Case/Qty	Part#
24-Well								
10.4mL		10-6	Yes	Yes	899410, 899405-	899410, 899405-1, 899403, 899406		931565-G-1X
10.4mL		10-6	Yes	Yes	Lid Included		20	931568
10.8mL		10-6	Yes	Yes	899410, 899405-	1, 899403, 899406	20	931569-G-1X
10.8mL		10-6	Yes	Yes	Lid Included		20	931571
96-Well								
500µL		non-sterile	*Yes	No	899410, 899403,	899406, 359747, 899405-1	50	9356045
650µL		non-sterile	Yes	No	899410, 899405-	1, 899403, 899406	50	931512B
2mL		non-sterile	*Yes	No	899410, 899405-	1, 899403, 899406, 359747	50	951657
2mL		10-6	*Yes	Yes	899410, 899405-	1, 899403, 899406, 359747	20	951657-S20
2mL		non-sterile	*Yes	No	899410, 899405-1, 899403, 899406, A210100		20	931130
2mL		10-6	*Yes	No	899410, 899405-1, 899403, 899406, A210100		20	931130-S
2mL		10-6	Yes	Yes	899410, 899405-1, 899403 899406, A210100		20	931133
2mL		10-6	Yes	Yes	Lid Included		20	931137
Filter P	lates							
Vol. Well	Well Shape	Sterility (SAL)	ANSI-SLAS	IW	Filter Membrane	Collection Plate	Case/Qty	Part#
24-Well Filt	er Plates							
10.8mL		non-sterile	Yes	No	25µm Polypropylene	931565-G-1X, 931568, 931569-G-1X, 931571	20	921550
~9mL		10-6	Yes	No	0.2µm Rapid Clear®	931565-G-1X, 931568, 931569-G-1X, 931571	20	921546
96-Well Filt	er Plates							
2mL		10-6	Yes	No	0.2µm Rapid Clear®	931130	20	921746
2mL		non-sterile	Yes	No	25µm Polypropylene	931130	25	931919
2mL		non-sterile	Yes	No	0.2µm PTFE	931130	20	921730
2mL		non-sterile	Yes	No	0.45µm PTFE	931130	20	921740
					0.2µm			

* Meets ANSI-SLAS plate dimensions	
------------------------------------	--

⁺ Irreversible Plate

Seals & Capmats

Description	Sterile	Plate Compatibility	Case/Qty	Part#
96-Well Capmat, For Wide Round Wells	No	951657, 9356045, 951657-S20, 931512B	50	359747
96-Well Capmat, For Square Wells	No	931130, 931130-S, 931133, 931137	100	A210100
Adhesive Foil Plate Seal	No	All Plates	100	899405-1
Pierceable Foil Heat Seal PCR compatible	No	All Plates	100	899403
Long-Term Storage Foil Heat Seal	No	All Plates	100	899406
Airporous Plate Seal For Growing Cultures	Yes	All Plates	100	899410
Well Plate Lid for use with 96- & 24-Well Plates	No	All Plates	100	981945
Well Plate Lid for use with 96- & 24-Well Plates	Yes	All Plates	100	981948

Vacuum Manifold

Sterile	ANSI-SLAS	Filter Membrane	Case/Qty	Part#
No	Yes	n/a	1	981802

