



## Introduction

The most critical aspects of reliable food contamination analysis are the reduction of interferences from the sample matrix and analyte recovery. Traditionally, SPE, SLE, Liquid-Liquid, syringe filtration, and centrifugation have been used to reduce matrix interference prior to LC/MS analysis. However, these techniques are time consuming, adversely impact recovery, require expensive consumables, and use large amounts of solvent (which need to be concentrated). Several studies were undertaken to investigate whether four different types of filter vial designs offered improved clean-up methods.

## Method I- Streamlined sample cleanup using combined dispersive solid-phase extraction in-vial filtration for analysis of pesticides and environmental pollutants in shrimp

L. Han, Y. Sapozhnikova, S.J. Lehotay, Anal. Chim. Acta (2014), http://dx.doi.org/10.1016/j.aca.2014.04.005

This process examines ways to enhance the overall method improvement of the analysis of pesticides and environmental pollutants in shrimp improvements for streamlining sample clean-up using dispersive and solid phase extraction in the Thomson eXtractor3D|FV will be compared traditional QuEChERS methodology. An effective way to reduce time and cost is to eliminate the centrifugation step and combine the SPE step filtration using the Thomson eXtractor3D|FV.

The sample matrix consisted of 42 diverse pesticides and 17 environmental contaminants in shrimp. Extracts were analyzed by both low press LC/MS/MS/MS.

## **Sample Preparation**

## **Extraction:**

- . 10g of spiked homogenized shrimp tissue, moisture content 86%. Atrazine-d5 was used as an Internal Standard. Samples were spiked at 10ng/g, 50ng/g and 100ng/g.
- 2. Add standards.
- 3. Vortex and allow to stand for at least 15 minutes.

### 4. Add 10mL of Acetonitrile to each sample.

- 5. Shake vigorously for 5 minutes.
- 6. Add 5g ammonium formate to each sample to induce phase separation.
- 7. Shake vigorously for 1 minute.
- 8. Centrifuge at room temperature for 2 minutes @ 4150rpm (3711rcf).

### Clean-up:

- 1.75mg of sorbent + 0.5mL of extract is added to eXtractor3D|FV shell.
- 2. Partially depress the Thomson eXtractor3D|FV p 0.2 µm PVDF membrane into the shell.
- 3. Shake for 30 seconds. 4. Completely depress the plunger into the shell to sample and analyze.

## htslabs.com htslabs.com/techcenter

## Method II- Multiclass multiresidue analysis of >100 veterinary drug residues in bovine tissues by filter-vial dispersive-SPE and LC-MS/MS

Steven J. Lehotay, Marilyn J. Schneider, and Alan R. Lightfield U. S. Department of Agriculture, Agricultural Research Service, Eastern Regional Research Center; 600 East Mermaid Lane; Wyndmoor, PA 19038; USA

High-throughput analysis is needed to meet demand for monitoring of veterinary drug residues in food animal tissues. The current veterinary drug residue monitoring method used by the USDA Food Safety and Inspection Service (FSIS) uses a combination of hexane-partitioning, dispersive-SPE, and solvent evaporation to achieve adequate cleanup for 20 mg equivalent sample injections to meet regulatory detection limit needs. This extra effort adds to the time and cost of the method and limits sample throughput. An improved method is developed and validated to streamline the sample preparation and LC-MS/MS method for identification and quantification of >120 veterinary drugs in bovine animal tissues for use in high throughput monitoring in the FSIS National Residue Program.



## **Experimental**

4. Inject 1µL in LC-MS/MS.

## Equipment

LC: Agilent 1100 MS/MS: AB Sciex 6500 Q-Trap, Electrospray ionization (positive/negative switching) Column: Phenomenex Kinetex C18 (50 x 3.0mm, 2.6μm) Mobile phase: (A) water; (B) MeCN both with 0.1% HCO<sub>2</sub>H Gradient: 2% A 100% B over 8.0min, hold for 2.7min Flow rate: 0.3 mL/min Column temperature: 40°C

## Validation and Results

- chemists

Method III- Routine Targeted Quantitation and Identification of Pesticide Residues using Triple Quadrupole LC-MS/MS and Advanced Scheduling of MRM **Transitions Detection of pesticides in filtered QuEChERS** extracts of Avocado Carrot, Grape, and Spinach

André Schreiber1, Julia Jasak2 1 AB SCIEX, Concord, ON (Canada); 2 AB SCIEX, Darmstadt (Germany)



# Improved Sample Clean-up Options for Contaminant Analysis for Vegetation, Meats & Seafood

Authors: Lisa Wanders, Sam Ellis, Joe Machamer, Thomson Instrument Company, Oceanside, CA

	Equipm	ent Co	nditions:	
	LPGC/MS/MS	– Agilent 7	390A	
and	• 220V fast oven h	neating upgra	de	
b	• Column – Restel film thickness	k non-coated	estrictor column & a Supelco SLBTM-5ms, 15m x 0.53mm x 1	١m
	• Vacuum outlet -	5.5m x 0.18m	m i.d.	
	Constant Flow –	2mL/min		
	• Carrier Gas – He			
	• Oven Temperatu	ure – 70°C for	.5 minutes	
o. Method I to the existing o with in-vial	• Injection – 5μL			
	LC/MS/MS – A	Agilent 1100	HPLC coupled to an Applied Biosystems API 3000	
sure GC/MS/MS and	MS/MS			
	<ul> <li>Electrospray ion</li> </ul>	ization in pos	itive mode	
	Source tempera	ture: was set t	o 525°C	
	• Column: Phenor	menex Revers	e Phase Prodigy ODS3 column, 150mm x 3.0mm x 5μm partic	le
	size			
the Thomson	Column Temper	ature: maintai	ned @ 30°C	
	• Flow rate: 0.3ml	_/min		
olunger with	• Mobile Phase:			
	A: 0.1% aqueous	formic acid		
o filter the	B: 100% Acetoni	trile		
	Time (min)	% <b>A</b>	% <b>B</b>	
	Initial	70	30	
	8.0	70	30	
	12.6	0	100	

## **Results:**

This approach to streamlining the QuEChERS protocol for the analysis of shellfish by combining the dispersive, sorbent and sample filtration into one vial, Thomson eXtractor3D|FV 0.2µm PVDF membrane, saves time, uses less solvent, and does not require special equipment. Table 1 shows the overall average recoveries of the 59 analytes in shrimp using different sorbents (n = 9 from triplicate spikes each at 10, 50, and 100 ng/g).Atrazine-d5 was used as the internal standard in both LPGC-MS/MS and HPLC-MS/MS. Recoveries of 13 of the 59 analytes were recovered at 100% while 42 of the 59 pesticides and contaminants tested were >70% with < 20% RSD independent of the sorbent used. The following pesticides were partially recovered depending on degradation, extraction partitioning factors, and specific sorbent used. Detection limits were < 5ng/g (with the exception of PCB's)

## Sample Preparation and Analysis

1. Weigh 2g homogenized tissue sample into 50mL tube. Add internal standard and spike solutions as needed. 2. Dispense 10 mL 4/1 (v/v) MeCN/water and shake 5min on a platform shaker, then centrifuge 5min at 3700rcf. 3. Transfer 0.4 mL extracts to filter-vial shell containing 25 mg C18, insert plunger halfway and shake 30s, then fully depress plunger to filter final extract into autosampler vial.

• Spikes made of 18 different tissue blank samples at 0X, 0.5X, 1X, and 2X levels (n=10 each) repeated 3 days by 3

• Matrix-matched and reagent-only calibration stds prepared at equiv. tissue levels of 0X, 0.25X, 0.5X, 1X, 2X, and 3X Internal standards were added, but not needed nor used • Method LOQs determined for spiked samples in matrix

## Updated method logistics compared to previous

FSIS Method Logistics (UPLC-TQD)	Ne
1 chemist was able to process 60 pre-homogenized samples in an 8-hr day for an overnight sequence	1 che
(longest step was 1 hr to evaporate MeCN)	(le st
No glassware to be cleaned afterwards	No g
Cost of materials $\approx$ \$3/sample (using bulk C18)	Wast
Waste = 10 mL hexane and 5 mL MeCN & two 50 mL, one 15 mL PP tubes, and AS vial	Revie took
Geis-Asteggiante et al., J. Chromatogr. A, 1258 (2012) 43-54 and Lehotay et al., Drug Test. Analysis, 4 (Suppl. 1) (2012) 75-90 and USDA-FSIS Analytical Chemistry Guidebook online	~~~



## **Experimental**

### Sample Prep:

1. SCIEX iDQuant<sup>™</sup> standards kit for Pesticide Analvsis

- 2. Store bought fruit & vegetables were extracted using European Standard Method 15662
- 3. Extracts were diluted 5x with water in Thomson Filter Vials, 0.45 µm **PVDF** membrane

### **Equipment Conditions:** • AB Sciex Triple Quad<sup>™</sup> 3500 with Turbo V<sup>™</sup> source and Electrospray Ionization

- Positive Polarit
- Column:
- Phenomenex Kinetex<sup>™</sup> Biphenyl 2.6 µm column

## Mobile Phase:

- Fast gradient of water/methanol with 5mM ammonium formate Flow rate: 0.5mL/min
- Injection 2 μ

## Result

- Average gain in sensitivity of 3x was observed
- Most pesticides had a detection limit of < 1ng/mL
- All pesticides had a detection limit of < 2ng/mL





system



both isomers of Dimethomorph. Examples of using the Flexible Window Width in a Scheduled MRM<sup>™</sup> Pro method: the window for Boscalid was set to 45 sec and Dimethomorph was detected using a wider window to detect both isomers together

Examples of using the Flexible Window Width in a Scheduled MRM<sup>™</sup> Pro method: the window for Boscalid was set to 45 sec and Dimethomorph was detected using a wider window to detect both isomers together



#	Analyte	MgSO, MgSO, + PSA		MaSO + PSA					
		+ filter	+ C18 + Z-Sep + filter	+ C18 + Z-Sep + CarbonX <sup>®</sup> + filter		#	Analyte		
		%Rec. (%RSD)							
1	Aldicarb*	99 (9)	104 (11)	103 (2)					
2	Aldicarb sulfone*	97 (6)	96 (9)	97 (6)	-	18	Ethoprophos*		
3	Atrazine	103 (3)	103 (4)	104 (4)		19	Fenthion		
3	Atrazine*	102 (4)	103 (4)	101 (6)		20	Fenthion sulfone		
4	Azoxystrobin	106 (9)	105 (14)	97 (12)		21	Imidacloprid*		
4	Azoxystrobin*	110 (4)	115 (3)	113 (3)	-	22	Lindane(γ-HCH)		
5	Carbaryl*	103 (4)	106 (4)	85 (7)		23	Linuron*		
6	Carbofuran	109 (12)	94 (17)	99 (12)	-	24	Methidathion		
6	Carbofuran*	104 (2)	109 (5)	109 (4)		25	Methidathion*		
7	Chlordane	108 (6)	102 (12)	104 (10)	-	26	Phosmet		
8	Chlorpropham	102 (11)	97 (14)	95 (14)		27	Phosmet*		
9	Chlorpyrifos	110 (7)	105 (11)	101 (11)	-	28	Pirimiphos-methyl		
10	Cyazofamid*	102 (5)	108 (5)	111 (4)		29	Pyriproxyfen		
11	Cypermethrin	113 (7)	107 (14)	102 (9)	-	30	Pyriproxyfen*		
12	o,p'-DDE	102 (7)	99 (11)	101 (11)		31	Tetrahydrophthalimide		
13	Deltamethrin	111 (10)	107 (15)	100 (16)	-	32	Tolclofos-methyl		
14	Diazinon	106 (9)	102 (12)	103 (11)		33	Triflumizole		
15	Dimethoate*	104 (6)	106 (7)	104 (7)	-	34	Triflumizole*		
16	Endosulfan sulfate	112 (8)	113 (6)	115 (8)		Enviro	nmental Contaminants		
17	Ethoprophos	103 (11)	91 (15)	91 (14)		35	PBDE 47	·	
			1		_	36	PBDE 99		



ew Method Logistics

nemist was able to process 60 pre-homogenized samples in 3 hours longest steps involved labeling tubes/vials, weighing, and preparing calibration standards)

glassware to be cleaned afterwards

ste = 10 mL MeCN and one 50 mL tube and an autosampler vial

iew of results for 135 drugs x 3 transitions x 67 injections (>27,000 data points)





## + • XIC from Data Scheduled 45.0, Gaussian resoluted Carbaryl ak on the state leke, dif Flufenacet State de la • XIC from Data Scheduled 88.1). Gaussian smoothed 2326 Methomyl • XIC from Data Scheduled 75.1). Gaussian smoothet Thiabendazole the second second

Fig. 2 Spinach sample showing 14.9 µg/kg Boscalid and 53.7 µg/kg

>3000 µg/kg Thiabendazole in orange Quantifier MRM above threshold Qualifier MRM triggered MRM above threshold Dynamic window extension 1. extension 2. extension Thiabendazole in orange (zoom) Quantifier MRM above threshold Qualifier MRM triggered MRM above threshold Dynamic window detection 1. extension 2. extension 

Fig. 3 Orange sample shows an example of MRM-triggered MRM and Dynamic Window Extension: the qualier MRM transition is automatically triggered when the quantier MRM transitions exceeds the threshold set in the Scheduled MRM<sup>™</sup> Pro method, the detection window is automatically extended if the MRM signal is above the threshold at the end of the detection window

Examples of MRM-triggered MRM and Dynamic Window Extension: the qualifier MRM transition is automatically triggered when the quantifier MRM transitions exceeds the threshold set in the Scheduled MRM<sup>™</sup> Pro method, the detection window is automatically extended if the MRM signal is above the threshold at the end of the detection window

## **Among-Day Reproducibility of Recoveries**

HorwitzRatio = HORRAT =  $RSD_{R}/(2C^{-0.1505})$ in which  $RSD_{R}$  is reproducibility and C is concentration (g/g)





Fig 4. Detection of pesticides in filtered QuEChERS extracts of avocado (A), carrot (C), grapes (G), and spinach (S) Quantitation and identification based on MRM ratios in MultiQuant<sup>™</sup> software, the example shows the side-by-side peak review for Boscalid with positive findings in grapes and spinach samples

	1	
gSO <sub>4</sub>	MgSO <sub>4</sub> + PSA	MgSO <sub>4</sub> + PSA
ilter	+ C18 + Z-Sep	+ C18 + Z-Sep
	+ filter	+ CarbonX®
		+ mter
Rec.		
RSD)		
2 (5)	101 (4)	101 (8)
8 (4)	108 (6)	100 (5)
9 (7)	111 (8)	103 (5)
0 (9)	106 (7)	98 (6)
7 (7)	82 (15)	83 (13)
9 (7)	106 (6)	87 (9)
7 (3)	106 (6)	106 (6)
7 (4)	112 (4)	111 (3)
8 (4)	106 (7)	91 (5)
7 (5)	111 (6)	99 (7)
1 (7)	106 (11)	105 (11)
4 (6)	104 (6)	95 (3)
5 (10)	114 (4)	98 (4)
0 (15)	90 (12)	93 (11)
1 (10)	98 (11)	90 (14)
4 (4)	95 (6)	90 (6)
7 (5)	103 (6)	99 (5)
3 (10)	102 (17)	93 (10)
3 (9)	94 (11)	70 (9)

#		Analyte	MgSO₄ + filter %Rec. (%RSD)	MgSO <sub>4</sub> + PSA + C18 + Z-Sep + filter	MgSO <sub>4</sub> + PSA + C18 + Z-Sep + CarbonX® + filter
37		PBDE 100	110 (10)	102 (8)	89 (8)
38		PCB 105	101 (8)	92 (15)	79 (7)
39		PCB114	99 (11)	98 (12)	82 (9)
40	,41	PCB118+123	98 (8)	91 (11)	80 (9)
42		PCB156	105 (8)	97 (14)	78 (14)
43		PCB157	100 (12)	96 (12)	76 (8)
44	-	PCB167	99 (11)	95 (12)	77 (8)
Int	ternal S	Standards			
45		Atrazine-d5	84 (5)	80 (10)	76 (6)
46	)	Atrazine-d5*	84 (9)	80 (5)	80 (6)
47		Fenthion-d6	86 (7)	81 (11)	72 (8)
* F	HPLC-M	S/MS results.			

Thomson Instrument Company is not affiliated with Phenomenx®, Ab Sciex®, USDA, Agilent® or any of their products..

Table 1: Overall average recoveries (and RSD) of the 59 analytes in shrimp using the 3D eXtractor filter vial d-SPE approach with different sorbents.

## Conclusions

The results clearly show the Thomson eXtractor3D|FV 0.2µm PVDF membrane approach to sample preparation using QuEChERS in an autosampler ready vial for partitioning, clean-up and filtration of shellfish is a fast and convenient method. This method lowers cost, solvent usage, and time. The Thomson eXtractor3D|FV yielded recoveries of 42 of the 59 pesticides and contaminants tested with >70% and < 20% RSD in shrimp. Future experiments will include optimization of dispersive and sorbent concentrations.

## Limits of Quantification (LOQs) for the Final Method



## **Conclusions**

• Our previous MRM for vet drug residues was streamlined and improved by using modern LC-MS/MS and eliminating hexane partitioning cleanup and solvent evaporation steps, and using filter-vial d-SPE for cleanup

• This method is more quantitative than our previous method

• LOQ <10 ng/g for nearly all drug analytes tested

• Sample throughput was 60 samples per day per 2 chemists for both sample prep and analysis/review of results

• Qualitative screening and identification results still need to be evaluated for method implementation by USDA-FSIS

Index	Sample Name	Sample Type	Dikter	Component Name	Masa jelo	Companient Group Name	Actual	Expected	Retention	Used	Ассанику	Calculated	las Ratie	Expected los Rute
7484	TH wondle	University	1.00	Filman 1	262.01 207.0	Benefit	11.4	10.68	1175	2	N/A	0.287	10000	13002
1455	TH minister	Life me	5.00	Email: 7	543.01 (40.0	Braceld	No.	10.98	HA.	12	RA	H=	NA	12201
116	70 (501	Unima	5.01	Administral 1	441.0 101	firmed.	Ne	10.08	84	-91	New.	14A	14.6	1.0000
1000	TH cant	(Jok-	东侧	Commod 1	343.07140.0	E-month	192	16.59	14.0	12	(ALA)	NiR	(g))	1.241
8275	Tri grapes	Unkresse	5.00	Boycarid 1	345.07.07.0	Boarakd	NA	10.55	10.59	121	No.2	11.344	3.0000	1.0000
8280	TH grapes	Linkraisen	5.00	Bosceld 2	343.07 140.0	Bascald	16.5	10.58	10.18	521	REA	17.718	0.2396	8:2400
8635	TH spitech	Unknown	5.00	.Bençaikd 7	345.0.) 307.9	Brazakd	100	10.58	10.60	90	NA	12.286	1.0000	1.0000
9882	TH spinach	Liekopan	6.00	Biocalid 2	343.0) \$85.0	Boostalid	Ach	10.50	10.59	191	NALA:	12,745	0.2540	0.2420
4123	Tri-sevals	1 Sol second	6.00	Discration I.	925/871	Brenahd.	18-16-	11.00	18.73	3	NEUL .	S'MI	1 A000	1.3054
3104	Tri Lonais	University	5.00	Bouwist?	363 2/ 120 0	Speakd	110-	10.55	169	12	NEW	184	1454	02430
德国	Tel armet	University	5.00	Expreshe 1	5450/387.6	Brend	14.12.	10.68	10.77	121	10.0	0/257	1,0000	3,0000
35.15	Di mor	Life man	A.08	Present 12.7	MAR 100 0	Benchil	11.4	0143	44	-141	N4.	114	NG	5,351
				- 11 - 2	The second second	-	TR. The second		- D2 T	( wanted		all'in TH		-
2 Boscar 120 200 100	id 1/ Standard Rod, Height 2 0 0 0 0	• 15 17 34Gd/J. (H av 59562 R. Area 50	1 evecado 01.ado - 8kis 2 97742 He 3037 - 2000 - 1300 - 1300 -	colid 1 (Unitrount. 17 ght 5 254a1 RT 4 1025	2 TH canol 4 convox - Bescaldo are fulli, Haright 1 3005 2000- 1002- 0	A 1 (Unit source) 1 UR, RT, MR min	21 TH preses TH prepers - Bosic Areas 1 530w4 He 3000 3000 1000		160 Comp	1 spinach mach Be 1 145e4 H 2000 - 1000 - 0	ecalid 1 (tieko weglet 2 25743 (0.60	25 Thomas (RT Anna KA Anna KA Anna KA	Esscalid T (L Height NGA, R7 200 000 000	tekroseni TerA mir

Fig 5. Quantitation and identification based on MRM ratios in MultiQuant<sup>™</sup> software, the example shows the side-by-side peak review for Boscalid with positive findings in grapes and spinach samples

### Conclusions

cleanup depending on the type of matrix.

• The developed method was applied to the quantitation and identification of pesticides in real food extracts.

• Following the European Standard Method 15662, different dispersive SPE kits from Phenomenex were used for sample

• Extracts were diluted 5 times with water to minimize possible matrix effects directly in the Thompson eXtreme|FV<sup>®</sup> 0.45 μm PVDF membrane part # 85541-500 and placed into the autosampler for LC-MS/MS analysis.