

Thomson 7L Flasks Grow 5L-6L:


Transient CHO/Stable CHO, HEK293, Insect Cells




2.5L (final volume)

5L (final volume)

Speed/Throw: 140 rpm with 25mm orbit
100 rpm with 50mm orbit



Scalable Protein Quality and Titer
in Transient CHO



6L, SF21 Insect Cells

Speed/Throw:
135 rpm with 25mm orbit

Thomson 7L: 50% volume+ more than other shake flasks, and can grow up to 30L/Shaker (eliminates bags)

Thomson 7L OG Flask Options



Optimum Growth® 7L Flask Scale
Consistently Across All Sizes
931117



Direct Fill to Bioreactors Optimum
Growth® 7L Multiport Flask
2' of ID 1/8" Transfer Tubing & Sampling
Port
931117-55-18



Direct Fill to Bioreactors Optimum
Growth® 7L Multiport Flask
2' of ID 1/4" Transfer Tubing & Sampling
Port
931117-55-14



Bidirectional Transfer Cap For Optimum
Growth® 7L Flask
931470-8

Thomson 7L: 50% volume+ more than other shake flasks, and can grow up to 30L/Shaker (eliminates bags)